

## Top

Welded corners are polished to a smooth, satin finish.

One-piece die-formed Type 304 stainless steel top provides a non-glare working surface.

Raised outer edge keeps fluids contained.

Sides are flanged down and hemmed to eliminate sharp edges.

Melamine-finished underside repels moisture for easy care and long life.

Convenient drawer opens from either side.

Complete frame supports top

## Base

Intersection of frame and cross-member is coped and welded.

Tubular U-shaped uprights are 1" OD, 16-gauge, Type 304, polished, stainless steel.

## Casters

Swivel casters, 5-in. dia. wheel; two with wheel brakes.

