

**CLASSIC PREP-PROCEDURE TABLES ARE DESIGNED WITH THESE QUALITY FEATURES:**

**Base**

Die-formed stainless steel, Type 304, wrap-around apron holds tub steady on frame.

Crossmember's ends are coped following the circumference of the tubing for strongest fit.

Junctions are fully circumference welded and polished to a smooth, crevice-free, satin finish.

Stainless steel, Type 304, satin-finished, 1-in. (2.54cm) round tubular crossmembers

Removable Prep-Procedure Rack has recessed patient positioning surface made of stainless steel round rod, welded at all intersections for strength. Has Electro-polish finish for easy cleaning .

Legs are fully circumference welded to the apron and polished to a smooth, crevice-free, satin finish.

Stainless steel, Type 304, satin-finished, 1.25-in. (3.17cm) round tubular legs

Adjustable floor levelers help keep table steady on uneven floors.

**Classic 6 Recessed-end Prep-Procedure Table,  
with 6 Single-Depth Tub and Removable Prep Rack  
Faucet and Drain sold separately, see Faucets Section.**

