

ALL SSCI LAMINATED ARISTOCRAT CAGES ARE DESIGNED WITH THESE QUALITY FEATURES:

Door

Stainless steel door spokes are deburred smooth to eliminate sharp ends.

Stainless steel, electropolished removable ambidextrous door may be remounted with hinges on the right side with no additional parts.

Door extends below cage floor to help keep paper liners in place.

Hinges

Heavy 12-gauge stainless steel hinges allows the door to operate smoothly.

Hinges are thru-bolted with a threaded connector to sides of cage using a 'T'-nut

Door fits snug to the cage body to eliminate areas where an animal's paws could be pinched.

.1875-in. (.48cm) diameter stainless steel vertical spokes are spaced 1-in. apart and welded at every intersection to crossmembers and frame.

.3125-in. (.79cm) diameter stainless steel door frame and cross-members are welded together for strength.

Latch

Stainless Steel Latch easily disassembles for sanitizing and adjusting.

.3125-in. (.79cm) diameter stainless steel latch keeper latches door securely.

Heavy 12-gauge stainless steel latch brackets.

Latch connector bar is notched for easy opening.

Heavy 12-gauge stainless steel self-latching door latch connector is designed to accommodate a padlock for greater security.

Body

All holes are precision-machined with a computer-controlled boring machine. Components are precision-milled in a computer-controlled milling machine.

Floor seams are sealed with an epoxy sealant. The sealant extends up the side joints for a few inches.

Distinctive floor design with fluid-retaining lip is standard on all Aristocrat Cages.

Plastic laminate is bonded to interior and front surfaces before assembly with a high-pressure press to create a superior moisture barrier.

Constructed of specially designed water-resistant .625 in. thick substrate Glued and mechanically-fastened at all joints for extra strength. All joints are bonded with a glue that is impervious to liquid.

