

# ALL OF OUR REGAL CASE CARTS ARE DESIGNED WITH THESE QUALITY FEATURES:

## Door

Removable double-panel stainless steel doors are designed for easy care: sloped top edges allow water to run off; bottom weep holes ensure effective drainage.

Doors open a full 270 degrees for loading and unloading Cart and latches in open position.

Door lifts off it's hinges for easy removal.

## Hinge

.3125-in. (.79cm) diameter stainless steel hinge pins.

Heavy 12-gauge stainless steel die-formed barrel hinges operate smoothly.

## Bumper

Full wrap-around vinyl bumper extends beyond the Cart sides and handle to protect facility doorways and walls.

## Body

One piece stainless steel top provides a smooth, crevice-free working surface.

Top surface may be used as an instrument table during procedures.

External handle extends the full width of the cart for easier maneuvering.

Integral door jamb on top, bottom and sides of door opening helps seal out environmental contamination.

Cart sides and back are made of a single piece of stainless steel formed with front flanges to create a vertical recessed door jamb.

Cart bottom is spot welded to sides and back for a smooth, easy to clean interior surface. Forward pitch and drain holes in the front corners ensure water can run out when cleaning.

## Frame

Four copper-nickel plated swivel casters with 5-in. (12.7cm) dia. Poly-Loc-Poly wheels have grease fittings and seals for long-lasting service. Two casters have swivel locks.

One-piece, welded, die-formed frame constructed entirely of stainless steel .

