

Vinyl gasket cushions back panel, providing a tight seal.

Full-perimeter stainless steel frame attaches rear full-view panel to cage for greater durability and a more finished appearance.

Features the same high quality as our standard Regal Cages.

Clear back panel available in polycarbonate or tempered glass.

Full-width and height clear back panel offers an unobstructed view allowing observation without direct interaction in post-op, boarding and adoption areas.

Vertical Molding Strips cover the space between two cages to keep liquids from seeping into hard-to-clean areas.

Corner Plates to join multiple cages together in the back to give a finished look.

