

CLASSIC 15 BATHING TABLE WITH TUBULAR STAINLESS STEEL BASE, RIGHTSIDE DOOR AND ROTATING WALK-UP RAMP

NEW!



**Never lift
a heavy, wet
animal in or out
of the tub
AGAIN!**

Efficient bathing equipment is an essential component of virtually every pet care facility. Whether your business specializes in regular preventative care, emergency services, boarding or grooming, you'll build customer satisfaction when you send your clients home from your facility bathed, groomed and sweet smelling.

Our Classic series Bathing Tables are designed to allow the most comfortable working height and allow for easy entry and exiting of large animals utilizing our rotating walk-up ramp in conjunction with our hinged self-latching access door. These Bathing Tables are available with ramp and door locations on either the left or right side. You'll never have to lift a heavy, wet animal in or out of the tub again. Imagine the convenience!

The entire Table — Tub, Splash Shields and Base — is made of Type 304 stainless steel, heli-arc welded for superior strength. All welds are ground smooth, and the entire tub assembly features a uniform, non-glare finish. One-piece Table design, exclusively manufactured by Suburban Surgical, incorporates .750-in. covered inside corners and .750-in. radius bends on all Tub and Splash Shields.

 **SUBURBAN SURGICAL CO., INC.**

For Product Information Call TOLL-FREE 1-800-323-7366 • (847) 537-9320 • Fax (847) 537-9061

Hours: 8:30am - 5:00pm CST Monday - Friday

Visit our website at: www.suburbansurgical.com

Bathing Table

One-piece Tub/Side Splash Shields with .750-in. (1.90cm) radius corners, and bottom lip that overlaps the tub edge keeps water away from walls.

Single lever faucet/mixing valve allows you to preset your maximum temperature setting.

Classic 15 Bathing Tables measure 60-in.W, 24-in.D, and 63-in.H overall.

.750-in. (1.90cm) radius bends on all sides and bottom give tub a smooth, easy-to-clean interior surface.

.750-in. (1.90cm) covered corners eliminate accumulation of dirt and debris.

Seams are fully crevice-free welded.

All welds are ground smooth, and the entire tub assembly features a uniform, non-glare finish.

3.5-in. (8.89cm) dia. opening for standard drain.

Bottom slopes 1 in. to drain.

Tub is 16-gauge, Type 304 Stainless Steel

Wall-mounted faucets and drains are sold separately, see Faucet Section.

Exclusive one-piece design incorporates .750-in. covered inside corners and .750-in. radius bends on tub and splash shields.



Base

Tubular legs and crossmembers are Stainless Steel, Type 304, satin-finished.

Crossmember ends are coped following the circumference of the tubing at each joint for strongest fit and are circumference welded.

Stainless Steel adjustable floor levelers help keep table steady on uneven floors.

All welds are ground and polished to a smooth, crevice-free, satin finish.

Access Door & Ramp

Rotating walk-up ramp has non-slip dura-grid surface which provides traction wet or dry.

Stainless steel hinged door swings in to allow for easy entry and exit. Includes heavy-duty, spring-loaded, self-latching positive door catch designed and used exclusively by Suburban Surgical.

Access door and rotating ramp are available on left or right side.

P/N	DESCRIPTION
13022-BH-JFKDDR	Classic 15 Bathing Table with Tubular Stainless Steel base, Right-side door and Rotating Walk-up Ramp
13023-BH-JFKDDR	Classic 15 Bathing Table with Tubular Stainless Steel base, Left-side door and Rotating Walk-up Ramp
109920	Wall-mounted, Single-lever with Riser, Flexible Extension hose and Straight Sprayer
024647-Z-Z-ZZ-ZZZZ-Z	Basket Strainer Drain Cover

SUBURBAN SURGICAL CO, INC.

For Product Information

Call TOLL-FREE 1-800-323-7366 • (847) 537-9320 •

Fax (847) 537-9061

Hours: 8:30am - 5:00pm CST Monday - Friday

Visit our website at: www.suburbansurgical.com