

## Cage Body

Convenient access ports may be used for IV lines or electronic monitoring equipment.

Plastisol-coated Cage Floors are included for patient comfort.

Removable wire grill and clear acrylic panel (included) slide into rear channel to create two compartments.

Built-in floor heaters provide even, steady heat to help maintain the patient's body temperature.

## Controls

Thermometer/hygrometer continuously measures the temperature and humidity inside the compartment.


Alarm sounds when oxygen falls below your desired concentration.

Oxygen meter monitors oxygen concentration.

Monitoring window lets you visually check when the carbon-dioxide filter needs to be changed.

Oxygen flowmeter regulates the flow of oxygen delivery.

Separate lights for each function indicate when function is on and operating.





Sturdy one-piece  
welded tubular  
frame.

Heli-arc welded  
joints seal out  
fluids, debris  
and bacteria.

All four Swivel  
casters with  
5-in. (12.70cm)  
diameter wheels  
have brakes for  
safety.

## Mobile Base for Intensive Care Unit



**Imperial/Regal Modular Cabinetry Base for Intensive Care Unit**