

Upper Grill

Stainless steel vertical spokes and crossmembers are welded together at every intersection.

Vertical spokes are made of .2375-in. (.603cm) diameter stainless steel.

Crossmembers are made of .312-in. (.792cm) diameter stainless steel.

Grill is welded to the tubular frame at all contact points to eliminate shifting, racking and rattling.

Frame

All corners have weldnut welded to the inside of the tubing to allow you to disassemble and reassemble pens without loss of structural integrity.

Frame is made of 18-gauge-1.25-in. (3.17cm) square Type 304 stainless steel tubing.

Panel, frame and grill feature a uniform, non-glare, easy-to-clean finish.

46.5-in. (118.11cm) high lower section is securely attached to the tubular frame.

Full-depth Floor-support Brackets hold floors securely in place.

Mitered corners and butted joints are heli-arc welded to seal out fluids and debris.

Leg and frame are made of solid, one-piece Type 304 stainless steel tubing.

Lower Panel

Lower section of End/Partition Panels is available in four styles:

- Single-Wall Stainless Steel
- Double-Wall Stainless Steel
- Fiberglass Reinforced Plastic in a choice of colors
- Water-resistant Laminate in a choice of colors

Adjustable Floor Levelers

Rear leveler is offset 14-in. (35.56cm) to accommodate a floor trough system.

Levelers adjust to fit irregularities in your floor.

