SUBURBAN SURGICAL COMPANY, INC. Stainless Steel Healthcare Products

800.323.7366 | 847.537.9320

275 Twelfth Street • Wheeling, Illinois 60090 • USA www.suburbansurgical.com

THE KEY WORD IS QUALITY.

SUBURBAN SURGICAL COMPANY specializes in designing, manufacturing and selling capital equipment to hospital professionals. Capital equipment purchases are not made often, therefore it is in your interest to be as knowledgeable as possible and purchase from a vendor who will help you at each step of the process.

For more than 60 years and three generations, we have designed and manufactured the highest quality equipment available for the health care market. Our goal is to provide our customers with the highest quality products, backed by the highest quality service. Our manufacturing capacity, warehouse capacity and support staff allows us to provide you with the quality equipment you need, when you want it, supported by the service you deserve.

Whether you're outfitting a brand-new facility, remodeling, moving, or simply looking to add some new equipment to enhance your services, SSCI is dedicated to meeting your needs. Call the SSCI sales department or your local authorized representative today and put us to work for you.

QUALITY PRODUCTS — OUR SPECIALTY

At SSCI, we design, engineer and manufacture quality health care equipment. We research the needs of the marketplace, design the best products to meet the needs, and manufacture all our products to exacting specifications to ensure they will provide excellent performance over the life of each item.

All that, plus the following distinctive features add up to make SUBURBAN SURGICAL the choice of medical care specialists nationwide.

STATE OF THE ART MANUFACTURING

We've made the investment in the best manufacturing equipment available, and will continue to do so as our company expands. Not only that, we have made a commitment to our manufacturing personnel to keep them up to date with the latest assembly processes and skills. This combination of quality machinery and experienced manufacturing staff makes SSCI the finest line of medical care equipment you can buy.

CONSTRUCTION DETAILS YOU WON'T FIND ANYWHERE ELSE

We listen to our customers to learn what's most important in the quality of the equipment we design, manufacture and sell. Take our extra-strong hinges on our Case Carts, or our 200+ pound capacity for our roll-out shelves. This is just an example of our exclusive design features that make our quality equipment the leader in the marketplace.

We want you to be happy with each and every item you purchase from us. That's why we go out of our way to provide the service you need before, during and after the sale.

SALES STAFF

Our experienced Sales Staff and Factory Authorized Representatives understand the needs of health care professionals, and can help you determine the best equipment for your use. Whether you buy from us direct or from a factory representative, you can call on our sales support staff to give you expert advice and answer any questions you may have about any of our products.

CUSTOM DESIGN SERVICES

SUBURBAN SURGICAL's Design Service is offered to customers after discussing and agreeing on a viable project with an SSCI Representative. After a methodical review of the customer's needs for their facility, the SSCI Representative will present the project to our staff of designers to create a set of approval drawings. The drawings ensure that SSCI will build exactly what the customer has discussed with the SSCI Representative. Once the customer has committed to purchase from SSCI, the Final Approval Drawings will be created and presented to the customer for their approval. Other Customer Requirements may be required before manufacturing time can begin.



IT'S OUR SPECIALTY — OUR DRIVING FORCE.

CUSTOMER SERVICE

Shipping Information

Once your order has been scheduled, our Customer Service staff will be available to answer any questions you may have including shipping information. Please review our Terms and Conditions, Customer Requirements and our Freight Receiving Procedures in the Business Procedure Section.

Installation Support

All our products are designed for quick and easy set up and installation. Detailed drawings are included with your products upon delivery. Our Customer Service staff can also answer any questions you may have during the installation process.

• Service after the Sale

Our Customer Service staff will be happy to help you with any questions or problems that arise after your equipment has been installed in your facility. Our goal is to be a long-term partner with your business. We're always just a phone call away.

WARRANTY

SUBURBAN SURGICAL COMPANY, INC. warrants the original purchaser that all equipment manufactured by Suburban Surgical Company, Inc. will be of the highest standards in material and workmanship. All equipment manufactured by Suburban Surgical Company, Inc. will be warrantied for a period of one (1) year from the date of shipment from the factory.

Components and Casters Suburban Surgical Company, Inc. purchases from other manufacturers will be covered by the respective manufacturers' Warranty.

Warranties will not apply if it is determined by Suburban Surgical Company, Inc. that the equipment became defective due to an accident, misuse, abuse, or alteration. Warranties do not include freight charges for replacement or repair.

WARRANTY AND NON-WARRANTY REPAIR POLICY:

• Please call our Customer Service Department at 1-800-323-7366, ext. 3518. 8:30AM - 5:00PM CST Monday - Friday.

Step 1: Discuss what is wrong with the product.

Step 2: SSCI will issue a Return Merchandise Authorization Number. The Return Merchandise Authorization Number must appear on the outside of shipping carton.

Step 3: Customer must package products in original shipping cartons and materials or in other appropriate shipping materials to prevent the products from being damaged in transit. If the products are damaged in transit and SSCI determines that the products were packaged improperly, SSCI will file a freight claim and call the customer with an estimate for repair. The customer can wait for the outcome of the freight claim or authorize the repair of the freight damages.

- SSCI will call the customer and advise if the repair is covered by the warranty. If the repair is not covered by warranty, SSCI will provide an estimate for the repair and ask for authorization to proceed.
- SSCI does not have products to loan to customers while their product is being repaired. By following the above procedure, SSCI will be able to make any necessary repairs in the minimum amount of time possible.

Warranty for Merchandise Manufactured of Stainless Steel is void if it is determined if they where subjected to harmful cleaning agents. See 'Guide to Care and Cleaning of Stainless Steel', packaged with each new Stainless Steel Product.

CLOSED CASE CARTS

Case Carts are ideal for keeping all the supplies needed for individual procedures close at hand. Simply load it with sterile instruments, surgical drapes, sutures and other supplies you need and wheel it to the surgical suite. When the procedure is finished, used items can be placed inside for return to facility's dirty area for disposing of waste and forwarding of reusable items to facility's sanitizing area where instrument's are cleaned and autoclaved. The Case Cart is restocked with sterile instrument packs and items needed for the next surgical procedure or emergency.

The quality of SUBURBAN SURGICAL's Stainless Steel Closed Case Carts quite simply can't be beat. We start with heavy-gauge Type 304 stainless steel, formed, heli-arc welded and polished to a smooth, even finish. We then incorporate a host of construction features shown above designed for industrial-strength durability and exceptionally easy care.

SINGLE DOOR CLOSED CASE CARTS

Single door is hinged left. Full-width external handle on right side. Cart can be configured with one or two optional shelves.

 105582-00
 77.79" W x 138.11" H x 69.53" D

 105587-00
 77.79" W x 100.01" H x 69.53" D

DOUBLE DOOR "OVERLAPPING" CLOSED CASE CARTS

 105602-00
 41.750" L x 27.375" W x 39.375" H

 105618-00
 53.750" L x 27.375" W x 39.375" H

STATIONARY & PULL-OUT WIRE SHELVES

Heavy-duty Type 304 stainless steel shelves feature continuous round rod construction welded at all intersections for extra strength. Turned-up back edge keeps items from falling off back when shelf is rolled out, and turned-down front edge eliminates a front lip to make it easier to load and unload items. The Pull-Out Shelf suspension tracks fit easily into shelf brackets inside cart. Shelf rolls out to 14.75-in. past Cart face, and positive stop prevents accidental removal. Height of shelves can be adjusted easily.

	Stationary Wire Shelf	Adjustable Roll-Out WIRE Shelf	Adjustable Roll-Out SOLID Shelf	Adjustable Roll-Out Perforated Shelf
Single Door Regal Case Cart	220300-00	220330-00	220340-00	220350-00
Double Door Regal Case Cart with Overlap	220300 00	220330 00	220310 00	220330 00
	220303-00	220331-00 220333-00 220335-00 220337-00	220341-00 220343-00 220345-00 220347-00	220351-00 220353-00 220355-00 220357-00
Double Door Regal Case Cart with center post				
care with center post	220307-00	220332-00 220334-00 220336-00	220342-00 220344-00 220346-00	220352-00 220354-00 220356-00

CARD HOLDERS AND ENGRAVED PLATES

Card Holders and Engraved Plates must be specified when ordering Case Carts.

600013

Open Top Edge, 5" x 5" x 3"



Stainless steel shelf runner has smoothrolling .750-in. (1.91cm) diameter polyethylene bearings.

> Retaining lip on sides and back keeps supplies from sliding off shelf when pulling shelf forward.





STAINLESS STEEL UTILITY/EQUIPMENT CART

Each Utility/Equipment Cart features a 36" L x 20" W Top that stands 34" H from the floor and one or two shelves. Both the Top and Shelves have a 3-sided Retaining Lip.

108110-00	Stainless Steel Utility/Equipment Cart with Shelf , 25.75" of open space between shelves, $39.809" \times 27.779" \times 40"$

108110-01 Stainless Steel Utility/Equipment Cart with 2 Shelves

STAINLESS STEEL TOTE-AL OPEN CASE CARTS

Our Open Case Carts are designed with a flat top and with or without a 3-sided Retaining Rail which prevents stored items from falling off the top while transporting. Total wrap around Bumper is mounted with brackets to provide clearance between the Bumper and Chassis providing adequate drainage. The Bumper prevents damage to facility walls.

105431-00	Overall Size without Retaining Rail 29.375" L x 27.875" W x 39.375" H
105433-00	Overall Size without Retaining Rail 35.375" L x 27.875" W x 39.375" H
105431-01	Overall Size with Retaining Rail 29.375" L x 27.875" W x 41" H
105433-01	Overall Size with Retaining Rail 35.375" L x 27.875" W x 41" H

INSTRUMENT TABLES



INSTRUMENT TABLES

Choose our stainless steel Instrument Tables to keep your instruments organized and right at hand during procedures. Use it to lay out instruments, then easily roll it to where instruments will be used. Constructed of Type 304 stainless steel for exceptionally easy care and long life. Choose from two styles; with a bottom cross brace or a shelf. Both styles are available in a choice of lengths and either with or without a drawer mounted beneath the top to offer greater storage capacity.

With Cross Brace, Overall Size

108410-00	20"L x 16"W x 34"H
108414-00	36"L x 24"W x 34"H
108415-00	48"L x 20"W x 34"H
108416-00	48"L x 24"W x 34"H
108417-00	60"L x 24"W x 34"H

Instrument Table With Shelf, Overall Size

108475-00	20"L x 16"W x 34"H
108479-00	36"L x 24"W x 34"H
108480-00	48"L x 20"W x 34"H
108481-00	48"L x 20"W x 34"H
108481-00	48"L x 24"W x 34"H
108482-00	60"L x 24"W x 34"H



INSTRUMENT TABLE WITH 3-SIDED RETAINING RAIL AND DRAWER

Keep your instruments organized and right at hand during procedures. Use it to lay out instruments, then easily roll it to where instruments will be used. Constructed of Type 304 stainless steel for exceptionally easy care and long life. Choose from two styles; with a bottom cross brace or a shelf. Both styles are available in a choice of lengths and either with or without a drawer mounted beneath the top to offer greater storage capacity.

108418-03 72.750"L x 24.375"W x 35.50"H Overall Size

I.V. AND MAYO STANDS

FOOT-OPERATED IV STANDS

Our I.V. Stands are designed for sturdiness and unfailing operation. The wide-spread, low center, heavy base design insures non-tipping stability. The Base Upright is heli-arc welded to the Base creating a one-piece, non-wobbling I.V. Stand. Smooth-rolling casters make for easy mobility. Hanger height is easily adjusted from 68"H to 108"H by simply raising to the desired height, automatically locking in place. Depressing the Foot-Operated release allows the Telescoping Tube to be lowered.

108180-00 2-Hook 108180-01 4-Hook 108180-02 6-Hook

HAND-OPERATED IV STANDS

Our I.V. Stands are designed for sturdiness and unfailing operation. The wide-spread, low center, heavy base design insures non-tipping stability. The Base Upright is heli-arc welded to the Base creating a one-piece, non-wobbling I.V. Stand. Smooth-rolling casters make for easy mobility. The knob-operated collet is designed to hold the Telescoping Tube at the desired Hanger height. Adjusts from 68"H to 108"H.

108200-00 2-Hook 108200-01 4-Hook 108200-02 6-Hook

FOOT-OPERATED MAYO STANDS

Our Mayo Stands are designed for sturdiness and unfailing operation. The wide-spread, low center, heavy base design insures non-tipping stability. The Base Upright is heli-arc welded to the Base creating a one-piece, non-wobbling Mayo Stand. Smooth-rolling casters make for easy mobility. Tray height is easily adjusted from 38"H to 61"H by simply raising to the desired height, automatically locking in place. Depressing the Foot-Operated release allows the Tray to be lowered.

Regular	Tray Size (WxD)	Overall Size (WxD)
108226-00	19.125" x 12.625"	19.50" x 17.50"
108228-00	21.250" x 16.250"	21.625" x 17.50"
Emergency	Tray Size (WxD)	Overall Size (WxD)
108231-00	19.125" x 12.625"	19.50" x 17.50"
108233-00	21.250" x 16.250"	21.625" x 17.50"



HAND-OPERATED MAYO STANDS

Knob-operated collet positively locks top at the desired height from 39-in. to 62-in. (99.06cm 157.48cm) high. Simply loosen knob to raise or lower top, then tighten to hold securely in place.

Regular	Tray Size (WxD)	Overall Size (WxD)
108246-00	19.125" x 12.625"	20.875" x 17.50"
108248-00	21.250" x 16.250"	23.000" x 17.50"
Emergency	Tray Size (WxD)	Overall Size (WxD)
108251-00	19.125" x 12.625"	20.875" x 17.50"
108253-00	21.250" x 16.250"	23.000" x 17.50"

Additional Mayo Trays

021905P	19.125" x 12.625"	48.58" x 32.07"	19.50" x 17.50"	49.53" x 44.45"
021906P	21.250" x 16.250"	53.98" x 41.28"	21.625" x 17.50"	54.93" x 44.45"

SOLUTION STANDS

MOBILE SOLUTION STANDS

Convenient Mobile Solution Stands are designed to keep 7-quart(s) (6.62-liter) round basin close at hand. Sturdy frame is made of 1-in. (2.54cm) diameter Type 304 stainless steel tubing, heli-arc welded at all joints for extra strength. Welds are buffed and polished to a bright shine. Comes complete with one or two 7-quart (6.62-liter) round basin(s) and four 3-in. (7.62cm) diameter non-marking, swivel casters.

108270-00	Stand with Cross Brace; one 7-qt. (6.62-liter) Round Basin, 13.88" W X 34.00" H X 13.88" D in.
108271-00	Stand with Shelf; one 7-qt. (6.62-liter) Round Basin 13.88" W X 34.00" H X 13.88" D in.
108272-00 29.88" W X	Stand with Cross Brace; two 7-qt. (6.62-liter) Round Basin < 34.00" H X 13.88" D in.
108273-00 29.88" W X	Stand with Shelf; two 7-qt. (6.62-liter) Round Basin < 34.00" H X 13.88" D in.
021758P Additional	7-quart (6.62-liter) Round Basin

RECEPTACLES & KICK BUCKETS

KICK BUCKET

An indispensable component of any surgical suite. Heavy-duty stainless steel bucket may be used under collection holes on operating tables or as an extra or replacement bucket for Solution Stands and Kick Bucket. Swinging handle folds out of the way.

021876P 13 quart (12.30-liter) Bucket

MOBILE SPONGE RECEPTACLE

Stainless steel base features a durable heli-arc welded frame with molded rubber pads and four 2-in. (5.08cm) diameter non-marking, swivel casters. Securely holds basin to keep it close to your work area. Vinyl wrap-around bumper protects frame and facility walls from damage. Available in a 7-quart (6.62-liter) round basin or a 10-quart (9.46-liter) oval basin.

108150-00	Base with 7-quart (6.62-liter) Round Basin
021758P	Additional 7-quart (6.62-liter) Round Basin
108145-00	Base with 10-quart (9.46-liter) Oval Basin
021714V	Additional 10-quart (9.46-liter) Oval Basin

MOBILE KICK BUCKET

Durable heli-arc welded stainless steel frame with molded rubber pads holds a 13-quart (12.30-liter) stainless steel bucket securely. Non-marking 2-in. (5.08cm) diameter casters swivel for superior mobility. Vinyl wrap-around bumper protects frame and your facility's furniture and walls from damage. Measures 13.50-in. diameter x 14-in. high (34.29cm diameter x 35.56cm high) overall.

108140-00 Includes 13-quart (12.30-liter) Bucket


OVER INSTRUMENT TABLES

OVER INSTRUMENT TABLE

Use this Table to keep surgical instruments, supplies, even monitoring equipment conveniently close at hand during procedures. Manufactured entirely of durable Type 304 stainless steel, our Over Instrument Table features a mobile telescoping base for versatility. Plus, our conveniently-designed crank operates from either side of the Table to raise top height from 39-in. to 62-in. (99.06cm to 157.48cm).

Table top surface measures a full 40-in. long x 23-in. wide (101.60cm long x 58.42cm wide) to hold an assortment of supplies. And the raised back rail keeps equipment from sliding off accidentally. Rolls smoothly on 3-in. (7.62cm) diameter non-marking casters to allow you to move Table into position easily. Two casters have brakes for safety.

108430-23 Top Size 36"L x 23"W, Height Adjusts 39" - 62", Clearance Between Legs 29"



33-in. (83.82cm) wide clearance between legs allows Table to easily straddle operating tables.

Caster brakes keep Table safely in position.

UTILITY TABLES

These Utility Tables feature built-in drawers to keep your instruments and treatment supplies clean, organized and easily accessible wherever you need them within your facility. Type 304 stainless steel construction offers long-lasting durability with easy care. Drawers are a full 4 in. (10.16cm) deep and feature a friction latch to stay closed while Table is being moved. Ball-bearing slides allow drawer to open fully, yet feature a positive stop to prevent drawers from being removed accidentally.

SINGLE & DOUBLE WIDE UTILITY TABLES

Our Single Wide Utility Tables are available in two sizes and are designed for long-lasting durable storage that is transportable. These sturdy tables can be used for a multitude of applications. Each Utility Table features a 3-sided Retaining Rail, 4" deep drawer(s) with friction latches and full extending ball bearing slides allowing drawer(s) to fully open to a positive stop preventing accidental removal. Most models also feature a Shelf.

108495-00 Single Wide Utility Table 20.750"L x 16.375"W x 35.656"H Overall Size Double Wide Utility Table with Shelf 36.750"L x 20.375"W x 35.656"H Overall Size 108520-00

DROP LEAF UTILITY TABLE

Drop Leaves (-03 models) are available for each Single Wide and Double Wide Utility Tables, adding 20" to the Top's length. The 3-sided Retaining Rail is removed when Drop Leaves are included with the Table.

Utility Table, One Drawer, with Drop Leaf 20.750"L x 16.375"W x 35.656"H Overall Width 40.750" 108495-03 108503-03 33.750"L x 18.375"W x 35.656"H Overall Width 53.750"

UTILITY STANDS

Our Utility Stands feature a 20"L x 16"W Top that stands 34"H from the floor and one Shelf. Both the Top and Shelf have a 4-sided Retaining Rail.

108300-00 SS Utility Stand, 20.750" L x 16.750" W x 36" H Overall SS Utility Stand with I.V. Post, I.V. Post adjusts from 68" H to 108" H, 20.750" L x 108301-00 16.750" W Overall

FOOT STOOLS

STAINLESS STEEL FOOT STOOLS

The Deck's non-slip surface and reinforced construction combine to create a stable, rigid platform. The "flared-out" unitized Leg design prevents tipping.

Single Step Fo	ot Stool	Overall Size
108330-00	16"L x 1	2"W x 8.25"H
108331-00	30"L x 1	2"W x 8.25"H

Two Step Foot S	Stool	Overall Size
108332-00	16"L x 2	4"W x 16.50"H
108333-00	30"L x 2	4"W x 16.50"H

FOOT STOOLS WITH HAND RAILS

SUBURBAN SURGICAL foot stools feature rugged sturdiness. The Deck's non-slip surface and reinforced construction combine to create a stable, rigid platform. The "flared-out" unitized Leg design prevents tipping. The welded-on Hand Rail creates added stability and security.

 Single Step Foot Stool, Overall Size

 108330-01
 16"L x 13.75"W x 38.625"H

 108331-01
 30"L x 13.75"W x 38.625"H

 Two Step Foot Stool, Overall Size

 108332-01
 16"L x 25"W x 47"H

 108333-01
 30"L x 25"W x 47"H

STACK AND LOCK FOOT STOOLS

Our Stack and Lock Foot Stools feature rugged sturdiness. The non-slip Vinyl Surface top, solid Stainless Steel construction, and Rubber Floor Tips combine to create a stable, rigid, transportable (and stackable) stool. Each Stack and Lock Foot Stool is designed to lock together with another Stack and Lock Foot Stool either side to side or front to back and be stacked on top of another. The handy Transport and Storage Cart holds up to 6 Stack and Lock Foot Stools.

Stack and Lock Foot Stool, Overall Size 108335-00 18"L x 12"W x 5"H

Transport & Storage Cart for Stack and Lock Foot Stool, Overall Size 222763-00 21.097"L x 12.150"W x 39"H

SCRUB SINKS

SURGEON SCRUB SINKS

Our 16-gauge, Type 304 Stainless Steel Surgeon Scrub Sinks are fully heli-arc welded into one piece to create the strongest Surgeon Scrub Sink available. Welds are continuous, crevice-free and ground smooth. Coved corners eliminate areas where dirt and debris accumulate. Radius bends on bottom and all sides and a uniform, non-glare finish results in a tub interior that is easy to clean. Scrub Sinks feature a 10.50-in. (26.67cm) deep basin, a 12-in. (30.48cm) high back splash. The Surgeon Scrub Sink is designed to accept a back splash-mounted goose neck faucet and 3.5-in. (8.89cm) drain. Stainless Steel Wall-Mounted Brackets are included.

A clean surgical scrub area is an essential part of sanitary, bacteria-free procedures. We make it easy for you to scrub up with these top-quality all stainless steel Surgeon Scrub Sinks.

Both sizes of Surgeon Scrub Sinks are designed to accept a gooseneck faucet. Faucets can be controlled by three types of hands-free valves; foot-operated, knee-operated and electronic-motion sensor.

13320-00-EPDMCX	Single-Station	30" W x 22.56" H x 19" D
13321-00-HJDMCX	Dual-Station	48" W x 22.56" H x 19" D

HANDS FREE FAUCETS

Rigid Gooseneck Spout is controlled by Knee-Action Valve. Adjustable valve mixes hot and cold water supply together. Bracket mounts valve on wall at the height you desire. Valve's body is chrome-plated brass. Purchase Spout and Valve together.

024642-Z-R-ZZ-ZZZZ-Z Rigid Gooseneck Spout with Aerator 024645-Z-Z-ZZZZZ-Z Knee-action Faucet Valve

ELECTRONIC SENSOR-ACTIVATED VALVE

Rigid Gooseneck Spout with Aerator is controlled by an Electronic Sensor Valve. Adjustable valve mixes hot and cold water supply together. Simply place your hands in front of the sensor to turn the water on. Powered by 120V/240V AC transformer that plugs into a wall outlet or DC requiring four "AA" cell batteries

024652-Z-R-ZZ-ZZZ-Z Rigid Gooseneck Spout and Valve

FOOT-ACTIVATED VALVE

Rigid Gooseneck Spout is controlled by Foot-Action Valve. Valve mounts to wall, positioning pedals just above floor. Hot and cold water supply are controlled by individual foot pedals. Valve's body is chrome-plated brass. Purchase Spout and Valve together.

024642-Z-R-ZZ-ZZZZ-Z Rigid Gooseneck Spout with Aerator 024646-Z-Z-ZZZZZZZZ Wall-mounted Foot-pedal Faucet Valve

BASKET STRAINER DRAIN COVER

Fits the drain in all SSCI tubs and sinks. Stainless steel construction stands up to years of use.

024647-Z-Z-ZZ-ZZ-Z Basket Strainer Drain Cover











SUBURBAN SURGICAL CO., INC.'S COMMITMENT TO THE RECOVERY OF THE UNITED STATES ECONOMY

Suburban Surgical Co., Inc. is proud to be a contributing business within the United States for more than 60 years. As a third generation family owned business, we have and continue to invest the skills and resources of our company, into our markets and economy. Not only are we well established in the animal laboratory housing, veterinary housing and operatory room markets, but also in the community through economic and charitable activities as well.

PRODUCT

Being the only company selling into the animal laboratory housing, veterinary housing and operatory room markets, we integrate in the value of our products: equipment, materials and components which directly support the jobs of both U.S. companies and U.S. citizens. Our primary supplier focus for materials used in the manufacture of our products is to purchase from U.S. based companies whenever possible.

- Total value of materials used to manufacture Suburban Surgical Co., Inc., products is typically in excess of 85% from U.S. based companies
- Tons of Stainless Steel purchased primarily from U.S. based companies exceeds 750,000

OPERATIONS

As a U.S. based company, Suburban Surgical Co., Inc. generates jobs and activities of various businesses. These businesses are directly related to our products and the support of our business.

Jobs, Operation Related: Marketing, sales, managing, design/engineer, manufacture, packaging, shipping and customer support.

Jobs, Product Related: Suppliers for raw materials, purchases of component parts, transportation of products are examples of support from U.S. based companies primarily.

SUBURBAN SURGICAL CO., INC, FACTS:

- Suburban Surgical Co., Inc. was established in the 1940's
- Suburban Surgical Co., Inc. was incorporated in 1971
- Suburban Surgical Co., Inc. is based in Wheeling, Illinois 60090
- Suburban Surgical Co., Inc. purchases from over 100 U.S. based companies
- Suburban Surgical Co., Inc. employs over 140 full time U.S. citizens

LOGISTICS

Suburban Surgical Co., Inc. is committed to the design, manufacture and timeliness of delivery to all of our customers.

- Operational facilities exceeding 200,000 square feet
- On time deliveries exceeding 95% on time or earlier
- Significant inventory of the most common products and components
- Inventory management systems which consistently monitor and replenish inventory

INDUSTRY AND COMMUNITY CONTRIBUTIONS

- Financial and In-Kind support of all local, regional and national industry organizations
- Direct employee involvement both in country, community and industry
- Financial charitable contributions

Suburban Surgical Co., Inc. is a U.S. based company that places an emphasis on supporting our country and growing the economy of the Unites States of America. Suburban Surgical Co., Inc. makes every effort to do so by employing American workers and facilitating American companies.

GREEN SUSTAINABLE DESIGN

The Economic Recovery Act specifically mentions "... grantees are encouraged to implement the following listed primary elements of sustainable design in federally funded facilities. Sustainability is the outcome of an integrated process of facility development and operation incorporating a balance of life-cycle costs, environmental impact, occupant health, safety, security and productivity."

Suburban Surgical Co., Inc. products are designed and manufactured with great attention to sustainability. Specifically, all of our products are designed and manufactured to the most stringent standards to ensure longer longevity. These standards not only reduce costs by extending our products "life-cycle" but also preserve our natural resources by reducing the need to repair or replace products on a more regular basis.



Main: 847.537.9320 Toll Free: 800.323.7366 Fax: 847.537.9061

www.suburbansurgical.com

Made in the USA