

“LEADER IN HOSPITAL CARE PRODUCTS FOR OVER 60 YEARS”



**SSCI**

**SUBURBAN SURGICAL Co., Inc.**

275 Twelfth Street, Wheeling, IL 60090

All intellectual material, including the text, graphics and other components of this catalog are Copyrighted by SSCI.  
Any use or duplication of these intellectual materials, without written consent by SSCI, is strictly forbidden.

# BASSINETS

## CLASSIC MATERIALS & CONSTRUCTION

### Frame

**Frame Members:**  
16 gauge, 1" dia. tubular stainless steel, type 304, bright buffed finish

**Basket Supports:** .375" dia. rod stainless steel, type 304, bright buffed finish

**Cross Members:** 16 gauge, 1" diameter tubular stainless steel, Type 304, bright buffed finish

**Cross Brace:** 16 gauge, 1" diameter tubular stainless steel, Type 304, bright buffed finish

**Shelf:** 18 gauge, stainless steel, Type 304, #4 finish

### Drawer

**Suspension Rails:** 16 gauge, stainless steel, type 304, #4 finish

**Housing:** Top: 20 gauge, stainless steel, type 304, #4 finish

**Back and Sides:** 20 gauge, stainless steel, type 304, #4 finish

**Spreader:** 20 gauge, stainless steel, type 304, 2B finish

**Slide Mounting Panel:** 20 gauge, stainless steel, type 304, 2B finish

**Drawer, Body and Front:** 22 gauge, stainless steel, type 304, #4 finish

**Pull:** dull chrome

**Slide Members, Housing/Drawer:** full extension, ball bearing, zinc plated

### Slide-Out Shelves

**Suspension Rails:** 16 gauge, stainless steel, type 304, #4 finish

**Retaining Rail:** .312" dia. rod stainless steel, type 304, bright buffed finish

**Shelf:** 18 gauge, stainless steel, type 304, #4 finish

**Legs:** "U"-shaped, 16 gauge, 1" dia. tubular stainless steel, type 304, bright buffed finish

**Casters:** 2" dia. non-marking wheel, swivel ball bearing, zinc plated

**Fasteners:** Machine Screws: stainless steel, type 18-8

**Lock Nuts:** stainless steel, type 18-8, nylon insert

**Rivets:** stainless steel, type 18-8

### Cabinet

**Suspension Rails:** 16 gauge, stainless steel, type 304, #4 finish

**Housing:** Top: 20 gauge, stainless steel, type 304, #4 finish

**Back and Sides:** 20 gauge, stainless steel, type 304, #4 finish

**Bottom:** 20 gauge, stainless steel, type 304, #4 finish

**Vertical Partition:** 20 gauge, stainless steel, type 304, #4 finish

**Removable Shelf:** 20 gauge, stainless steel, type 304, #4 finish

**Doors, Outer Panel:** 20 gauge, stainless steel, type 304, #4 finish

**Inner Panel:** 22 gauge, stainless steel, type 304, #4 finish

**Hinges:** stainless steel Pull: dull chrome

**Latch:** spring-loaded nylon



**Classic Bassinet with Cross Brace**

## Classic Bassinets are

designed for maneuverability, stability, functionality, and durability. Easy to clean Stainless Steel construction keeps risk of cross-infection at a minimum while maintaining a gleaming, high-polished appearance over time.

The optional, removable basket sits in either the flat or Trendelenburg position from either end. Bassinets are available with Drawers, Cabinets, Shelves, Slide-Out Shelves, or a combination of these, adding convenience and storage area.

All Bassinets are sold without the basket and the basket pad.



**Classic Bassinet with Cross Brace and 3 Compartment Storage Cabinet**



**Classic Bassinet with Cross Brace and 2 Door Cabinet**



**Classic Bassinet with Cross Brace and Side Opening Drawer**



**Classic Bassinet with Cross Brace and End Opening Drawer**

## Classic Bassinets with Cross Brace

P/N	Description	Overall Dimensions (LxWxH)			Cabinet/Drawer Dimensions (LxWxD)		
108000-00	Classic Bassinet with Cross Brace*	31.25"	15.75"	36.00"			
108000-01	Classic Bassinet with Cross Brace and 3 Compartment Storage Cabinet*	31.25"	15.75"	36.00"	30.00"	14.312"	19.00"
108000-02	Classic Bassinet with Cross Brace and 2 Door Cabinet*	31.25"	15.75"	36.00"	30.00"	14.312"	19.00"
108000-03	Classic Bassinet with Cross Brace and Side Opening Drawer*	31.25"	15.75"	36.00"	21.00"	12.750"	8.375"
108000-04	Classic Bassinet with Cross Brace and End Opening Drawer*	31.25"	15.75"	36.00"	26.00"	11.125"	8.375"

\* Basket not included.

# Classic Bassinets



**Classic Bassinet with Shelf**



**Classic Bassinet with Shelf and Side Opening Drawer**



**Classic Bassinet with Shelf and End Opening Drawer**

**Classic Bassinets** with the Slide-Out Shelf extends to a positive stop preventing accidental removal.



**Classic Bassinet with Cross Brace Slide-Out Shelf with 3-sided Retaining Rail**



**Classic Bassinet with Cross Brace Slide-Out Shelf with 3-sided Retaining Rail & Side Opening Drawer**

## Classic Bassinets with Shelf OR Cross Brace and Slide-Out Shelf

P/N	Description	Overall Dimensions (LxWxH)			Cabinet/Drawer Dimensions (LxWxD)		
108001-00	Classic Bassinet with Shelf*	31.25"	15.75"	36.00"			
108001-03	Classic Bassinet with Shelf Side Opening Drawer*	31.25"	15.75"	36.00"	21.00"	12.75"	8.375"
108001-04	Classic Bassinet with Shelf End Opening Drawer*	31.25"	15.75"	36.00"	26.00"	11.125"	8.375"
108002-00	Classic Bassinet with Cross Brace Slide-out Shelf with 3-sided Retaining Rail*	32.25"	15.75"	36.00"	62.25" slide out shelf extended		
108002-03	Classic Bassinet with Cross Brace Slide-out Shelf with 3-sided Retaining Rail and Side Opening Drawer	32.25"	15.75"	36.00"	21.00"	12.75"	8.375"
					62.25" slide out shelf extended		

\* Basket not included.



**Classic Stainless Steel Bassinet with Cross Brace  
Slide-Out Shelf with 3-sided Retaining Rail  
& 3 Compartment Storage Cabinet**



**Classic Stainless Steel Bassinet with Cross Brace  
Slide-Out Shelf with 3-sided Retaining Rail &  
2 Door Cabinet**



**Classic Bassinet with Shelf  
Slide-Out Shelf with 3-sided Retaining Rail**



**Classic Bassinet with Shelf  
Slide-Out Shelf with 3-sided Retaining Rail &  
Side Opening Drawer**

## Classic Bassinets with Shelf AND Slide-Out Shelf

P/N	Description	Overall Dimensions (LxWxH)			Cabinet/Drawer Dimensions (LxWxD)		
108002-01	Classic Bassinet with Cross Brace Slide-out Shelf with 3-sided Retaining Rail & 3 Compartment Storage Cabinet*	32.25"	15.75"	36.00"	62.25" slide out shelf extended 30.00" 14.312" 19.00"		
108002-02	Classic Bassinet with Cross Brace Slide-out Shelf with 3-sided Retaining Rail & 2 Door Cabinet*	32.25"	15.75"	36.00"	62.25" slide out shelf extended 30.00" 14.312" 19.00"		
108003-00	Classic Bassinet with Shelf Slide-out Shelf with 3-sided Retaining Rail*	32.25"	15.75"	36.00"	62.25" slide out shelf extended		
108003-03	Classic Bassinet with Shelf Slide-out Shelf with 3-sided Retaining Rail Side Opening Drawer	32.25"	15.75"	36.00"	62.25" slide out shelf extended 21.00" 12.75" 8.375"		

\* Basket not included.

# Imperial Bassinets

**Frame**

**Frame Members:** 16 gauge, 1" dia. tubular stainless steel, type 304, bright buffed finish

**Basket Supports:** .375" dia. rod stainless steel, type 304, bright buffed finish

**Cross Members:** 16 gauge, 1" diameter tubular stainless steel, Type 304, bright buffed finish

**Shelf:** 18 gauge, stainless steel, Type 304, #4 finish

**Casters:** 3" dia. non-marking wheel, swivel ball bearing, zinc plated

**Drawer**

**Suspension Rails:** 16 gauge, stainless steel, type 304, #4 finish

**Housing, Top:** 20 gauge, stainless steel, type 304, #4 finish

**Housing, Back and Sides:** 20 gauge, stainless steel, type 304, #4 finish

**Housing, Slide Mounting Panel:** 20 gauge, stainless steel, type 304, 2B finish

**Drawer, Pull:** dull chrome

**Slide Members, Housing/Drawer:** full extension, ball bearing, zinc plated

**Housing, Spreader:** 20 gauge, stainless steel, type 304, 2B finish

**Drawer, Body and Front:** 22 gauge, stainless steel, type 304, #4 finish

**Slide-Out Shelves**

**Suspension Rails:** 16 gauge, stainless steel, type 304, #4 finish

**Retaining Rail:** .312" dia. rod stainless steel, type 304, bright buffed finish

**Shelf:** 18 gauge, stainless steel, type 304, #4 finish

**Legs:** "U"-shaped, 16 gauge, 1" dia. tubular stainless steel, type 304, bright buffed finish

**Casters:** 2" dia. non-marking wheel, swivel ball bearing, zinc plated

**Cabinet**

**Suspension Rails:** 16 gauge, stainless steel, type 304, #4 finish

**Housing, Top:** 20 gauge, stainless steel, type 304, #4 finish

**Housing, Back and Sides:** 20 gauge, stainless steel, type 304, #4 finish

**Housing, Bottom:** 20 gauge, stainless steel, type 304, #4 finish

**Vertical Partition:** 20 gauge, stainless steel, type 304, #4 finish

**Removable Shelf:** 20 gauge, stainless steel, type 304, #4 finish

**Doors, Outer Panel:** 20 gauge, stainless steel, type 304, #4 finish

**Doors, Inner Panel:** 22 gauge, stainless steel, type 304, #4 finish

**Doors, Hinges:** stainless steel

**Doors, Pull:** dull chrome

**Doors, Latch:** spring-loaded nylon

## Imperial Stainless Steel Bassinets

Feature and extra-strong one-piece frame and casters set-in towards the center from the ends for increased mobility. Imperial Bassinets are available with Drawers, Cabinets, Shelves, Slide-Out Shelves, or a combination of these, adding convenience and storage area.

The Chassis is made from 16 gauge, 1" tubular stainless steel, type 304 and has 3" diameter casters with non-marking wheels with swivel. The Basket is .375" diameter rod stainless steel, type 304, and everything has a bright buffed finish.





**Imperial Bassinet**



**Imperial Bassinet with Side Opening Drawer**



**Imperial Bassinet with End Opening Drawer**



**Imperial Bassinet with 3 Compartment Storage Cabinet**



**Imperial Bassinet with 2 Door Cabinet**

## Imperial Bassinets

P/N	Description	Overall Dimensions (LxWxH)			Cabinet/Drawer Dimensions (LxWxD)		
108004-00	Imperial Bassinet*	31.25"	15.75"	36.00"			
108004-03	Imperial Bassinet Side Opening Drawer*	31.25"	15.75"	36.00"	21.00"	12.750"	8.375"
108004-04	Imperial Bassinet End Opening Drawer*	31.25"	15.75"	36.00"	26.00"	11.125"	8.375"
108004-01	Imperial Bassinet 3 Compartment Storage Cabinet*	31.25"	15.75"	36.00"	30.00"	14.312"	19.00"
108004-02	Imperial Bassinet 2 Door Cabinet	31.25"	15.75"	36.00"	30.00"	14.312"	19.00"

\* Basket not included.

# Imperial Bassinets



**Imperial Bassinet with Shelf**



**Imperial Bassinet with Shelf and Side Opening Drawer**



**Imperial Bassinet with Shelf and End Opening Drawer**

**Imperial Bassinets** with the Slide-Out Shelf extends to a positive stop preventing accidental removal.



**Imperial Bassinet Slide-Out Shelf with 3-sided Retaining Rail**



**Imperial Bassinet, Slide-Out Shelf with 3-Sided Retaining Rail & Side Opening Drawer**

## Imperial Bassinets with Shelf OR Slide-Out Shelf

P/N	Description	Overall Dimensions (LxWxH)			Cabinet/Drawer Dimensions (LxWxD)		
108005-00	Imperial Bassinet with Shelf*	31.25"	15.75"	36.00"			
108005-03	Imperial Bassinet with Shelf Side Opening Drawer*	31.25"	15.75"	36.00"	21.000"	12.750"	8.375"
108005-04	Imperial Bassinet with Shelf End Opening Drawer*	31.25"	15.75"	36.00"	26.000"	11.125"	8.375"
108006-00	Imperial Bassinet Slide-out Shelf with 3-sided Retaining Rail*	32.25"	15.75"	36.00"	62.250" slide out shelf extended		
108006-03	Imperial Bassinet Slide-out Shelf with 3-sided Retaining Rail Side Opening Drawer	32.25"	15.75"	36.00"	62.250" slide out shelf extended 21.000" 12.750" 8.375"		

\* Basket not included.





**Imperial Stainless Steel Bassinet with Slide-Out Shelf with 3-sided Retaining Rail & 3 Compartment Storage Cabinet**



**Imperial Stainless Steel Bassinet with Slide-Out Shelf with 3-sided Retaining Rail & 2 Door Cabinet**



**Imperial Bassinet with Shelf Slide-Out Shelf with 3-sided Retaining Rail**



**Imperial Bassinet with Shelf Slide-Out Shelf with 3-sided Retaining Rail & Side Opening Drawer**

## Imperial Bassinets with Shelf AND Slide-Out Shelf

P/N	Description	Overall Dimensions (LxWxH)	Cabinet/Drawer Dimensions (LxWxD)
1108006-01	Imperial Bassinet Slide-out Shelf with 3-sided Retaining Rail 3 Compartment Storage Cabinet*	32.25" 15.75" 36.00"	62.250" slide out shelf extended 30.000" 14.312" 19.000"
108006-02	Imperial Bassinet Slide-out Shelf with 3-sided Retaining Rail & 2 Door Cabinet*	32.25" 15.75" 36.00"	62.250" slide out shelf extended 30.000" 14.312" 19.000"
108007-00	Imperial Bassinet with Shelf Slide-out Shelf with 3-sided Retaining Rail*	32.25" 15.75" 36.00"	62.250" slide out shelf extended
108007-03	Imperial Bassinet with Shelf Slide-out Shelf with 3-sided Retaining Rail Side Opening Drawer	32.25" 15.75" 36.00"	62.250" slide out shelf extended 21.000" 12.750" 8.375"

\* Basket not included.

# Classic and Imperial Bassinet

## Materials

### Support Structure/Chassis, Classic style:

**Frame Members:** 16 gauge, 1" dia. tubular stainless steel, type 304, bright buffed finish

**Basket Supports:** .375" dia. rod stainless steel, type 304, bright buffed finish

**Cross Members:** 16 gauge, 1" dia. tubular stainless steel, type 304, bright buffed finish

**Cross Brace:** 16 gauge, 1" dia. tubular stainless steel, type 304, bright buffed finish

**Shelf:** 18 gauge, stainless steel, type 304, #4 finish

### Support Structure/Chassis, Imperial style:

**Frame Members:** 16 gauge, 1" dia. tubular stainless steel, type 304, bright buffed finish

**Basket Supports:** .375" dia. rod stainless steel, type 304, bright buffed finish

**Cross Members:** 16 gauge, 1" dia. tubular stainless steel, type 304, bright buffed finish

**Shelf:** 18 gauge, stainless steel, type 304, #4 finish

### Casters:

3" dia. non-marking wheel, swivel ball bearing, zinc plated

### Drawer Unit:

**Suspension Rails:** 16 gauge, stainless steel, type 304, #4 finish

#### Housing:

**Top:** 20 gauge, stainless steel, type 304, #4 finish

**Back and Sides:** 20 gauge, stainless steel, type 304, #4 finish

**Spreader:** 20 gauge, stainless steel, type 304, 2B finish

**Slide Mounting Panel:** 20 gauge, stainless steel, type 304, 2B finish

#### Drawer:

**Body and Front:** 22 gauge, stainless steel, type 304, #4 finish

**Pull:** dull chrome

**Slide Members, Housing/Drawer:** full extension, ball bearing, zinc plated

### Cabinet:

**Suspension Rails:** 16 gauge, stainless steel, type 304, #4 finish

#### Housing:

**Top:** 20 gauge, stainless steel, type 304, #4 finish

**Back and Sides:** 20 gauge, stainless steel, type 304, #4 finish

**Bottom:** 20 gauge, stainless steel, type 304, #4 finish

**Vertical Partition:** 20 gauge, stainless steel, type 304, #4 finish

**Removable Shelf:** 20 gauge, stainless steel, type 304, #4 finish

#### Doors:

**Outer Panel:** 20 gauge, stainless steel, type 304, #4 finish

**Inner Panel:** 22 gauge, stainless steel, type 304, #4 finish

**Hinges:** stainless steel

**Pull:** dull chrome

**Latch:** spring-loaded nylon

### Slide-out Shelves:

**Suspension Rails:** 16 gauge, stainless steel, type 304, #4 finish

**Retaining Rail:** .312" dia. rod stainless steel, type 304, bright buffed finish

**Shelf:** 18 gauge, stainless steel, type 304, #4 finish

**Legs:** "U"-shaped, 16 gauge, 1" dia. tubular stainless steel, type 304, bright buffed finish

**Casters:** 2" dia. non-marking wheel, swivel ball bearing, zinc plated

### Fasteners:

**Machine Screws:** stainless steel, type 18-8

**Lock Nuts:** stainless steel, type 18-8, nylon insert

**Rivets:** stainless steel, type 18-8

## Construction:

### Support Structure/Chassis, Classic or Imperial styles:

**Frame Members:** one piece "U"-shaped or "oblong"-shaped; bends are mandrill-bent for wrinkle-free bends

**Basket Supports:** formed to fully support bassinet basket and fully heli-arc circumference-welded to Cross Members; welds are buffed to a smooth, bright polished finish

**Cross Members:** ends are coped for uniform fitting joints and fully heli-arc circumference-welded to Frame Members; welds are buffed to a smooth, bright polished finish

**Cross Brace:** ends are coped for uniform fitting joints and fully heli-arc circumference-welded to Cross Member; welds are buffed to a smooth, bright polished finish

**Shelf:** one piece with front, back, and both sides hemmed; hemmed edges are flanged down creating a unitized top and double-thick angle frame; corners are notched to fit the leg contour; unitized top and double-thick angle frame's vertical flanges are fully heli-arc welded to the chassis' frame members (welds are concealed) (Note: In models where Shelf is used, the lower Cross Members and Cross Brace are eliminated.)

### Casters:

**Classic Style:** stem-mounted into caster socket with caster socket inserted into ends of Frame Member

**Imperial Style:** threaded-stem mounted through holes in Frame Member and secured with lock nut

### Drawer Unit:

**Suspension Rails:** formed angles; vertical flanges are fully heli-arc welded to chassis Frame Members; welds are buffed to a smooth, bright polished finish

#### Housing:

**Top:** one piece with back and both sides flanged down; front is flanged down and formed to create a recessed horizontal drawer jamb

**Back and Sides:** one piece wrapper; bottom edge is hemmed; top and bottom edges are flanged in; back and sides are formed with 90 degree bends; the front flanges of the sides are formed to create a vertical recessed drawer jamb; wrapper is inserted into the Top with back, sides, and front flanges spot-welded to Top's flanges

**Spreader:** formed channel; heli-arc welded to Sides' bottom flanges (welds are concealed)

**Slide Mounting Panel:** formed channel punched to accommodate housing slide members; heli-arc welded to Sides' top and bottom flanges (welds are concealed)

#### Drawer:

**Body:** one-piece punched and formed creating the drawer's inner front, back, sides, and bottom with radius formed bottom bends and coved corners; front, back, and both sides' top edges are flanged out; front, back and sides are spot welded

**Outer Front Panel:** one piece punched with top edge formed into a channel; bottom edge and side edges are formed into channels with flanges; Outer Front Panel's top edge channel overlaps Drawer's inner front panel top edge with bottom and side flanges spot welded to the Drawer Body's bottom and sides

**Pull:** fastened to Drawer Front with machine screws

**Slide Members, Housing/Drawer:** riveted to Housing's Slide Mounting Panel and Drawer's sides

**General:** Drawer Unit is inserted into Suspension Rails and heli-arc welded; welds are polished to a uniform finish

# Classic and Imperial Bassinet Accessories

## Cabinet:

**Suspension Rails:** formed angles; vertical flanges are fully heli-arc welded to chassis Frame Members; welds are buffed to a smooth, bright polished finish

### Housing:

**Top:** one piece with back and both sides flanged down; front is flanged down and formed to create a recessed horizontal door jamb

**Back and Sides:** one piece wrapper; back and sides are formed with 90 degree bends; the front flanges of the sides are formed to create a vertical recessed door jamb; wrapper is inserted into the Top with back, sides, and front flanges spot-welded to Top's flanges

**Bottom:** one piece with front, back, and both sides flanged down; Bottom is inserted into the wrapper with back, sides, and front flanges spot-welded to wrapper to create an edge-free cabinet interior

**Vertical Partition:** one piece with front, back, top, and bottom flanged; inserted into Housing and spot welded to Housing's top, back, and bottom

**Removable Shelf:** one piece; front, back, and both sides are hemmed and flanged;

### Doors:

**Outer Panel:** one piece formed with top, bottom, and both sides flanged

**Inner Panel:** one piece formed with top, bottom, and both sides flanged; inserted into outer panel; corners are fully heli-arc welded; welds are ground smooth and polished to a uniform finish; intermediate heli-arc welds along top, bottom, and both sides; welds are polished to a uniform finish

**Hinges:** continuous; heli-arc welded to Door's edge and Housing's door jamb; welds are polished to a uniform finish

**Pull and Latch:** fastened to Door with machine screws

**General:** Cabinet is inserted into Suspension Rails and heli-arc welded; welds are polished to a uniform finish

## Slide-out Shelves:

**Suspension Rails:** formed angles fitted with polyethylene slides and heli-arc welded to chassis Frame Members; welds are buffed to a smooth, bright polished finish

**Retaining Rail:** one piece, 3-sided, with formed support ends; intermediate supports are fully heli-arc circumference-welded to Rail; Rail is heli-arc circumference-welded to top of Slide-out Shelf and Legs; welds are buffed to a smooth, bright polished finish

**Shelf:** one piece with front, back, and both ends flanged down and formed into channels; front and back channels are fitted with polyethylene slides; 2 end corners are fully heli-arc welded; welds are ground smooth and polished to a uniform finish

**Legs:** bends are mandrill-bent for wrinkle-free bends; ends are capped and heli-arc circumference-welded; welds are ground smooth and polished to a uniform finish

**Casters:** threaded-stem mounted through holes in Frame Member and secured with lock nut

## Removable Acrylic Bassinet Basket and Pad

Removable Acrylic Bassinet Basket

### Removable Bassinet Basket Pad

Airfoam Pad Covered with Plastic Sheeting  
26"L x 12.50"W x 2" Thick



P/N	DESCRIPTION	OVERALL SIZE (LxWxH)		
220000-00	Acrylic Bassinet, No Card Holder	31.50"	16.00"	9.25"*
220000-01	Acrylic Bassinet, Card Holder at Head	31.50"	16.00"	9.25"*
220000-02	Acrylic Bassinet, Card Holder at Foot	31.50"	16.00"	9.25"*
220006-00	Removable Bassinet Basket Pad	26.00"	12.50"	2.00"

\* 8"H Foot End

For Product Information Call TOLL-FREE 1-800-323-7366 • (847) 537-9320 • Fax (847) 537-9061

Hours: 8:30am - 5:00pm CST Monday - Friday

Visit our website at: [www.suburbansurgical.com](http://www.suburbansurgical.com)