

***Stainless Steel Warming Cabinet***

“LEADER IN HOSPITAL CARE PRODUCTS FOR OVER 60 YEARS”



***SSCI***

***SUBURBAN SURGICAL Co., Inc.***

275 Twelfth Street, Wheeling, IL 60090

All intellectual material, including the text, graphics and other components of this catalog are Copyrighted by SSCI.  
Any use or duplication of these intellectual materials, without written consent by SSCI, is strictly forbidden.

**For Product Information Call TOLL-FREE 1-800-323-7366 • (847) 537-9320 • Fax (847) 537-9061**

**Hours: 8:30am - 5:00pm CST Monday - Friday**

Visit our website at: [www.suburbansurgical.com](http://www.suburbansurgical.com)

# WARMING CABINET

**PRECISION CRAFTED STAINLESS STEEL EQUIPMENT  
DESIGNED FOR LONG-LASTING PERFORMANCE**

## OUR STAINLESS STEEL WARMING CABINET INCLUDES THESE FEATURES:

### Cabinet

Double wall insulated  
Stainless steel, Type 304,  
satin-finished.

Stainless Steel adjustable  
wire shelf, allows positioning  
for blankets and/or  
solutions.

Cabinet is made to sit on a  
counter top or Mobile Cart.

Heating element and fan are  
enclosed in the bottom of  
the unit and covered with a  
perforated shelf.



### HEATING ELEMENT

Controller is also equipped with  
both a visual and audible  
alarm, to prevent overheating.

Temperature controller ranges  
from 90°-160°F. The displayed  
can be set to show either  
Fahrenheit or Celsius.

The keypad is conveniently  
located and features: On/Off,  
Set Temperature, Alarm Light,  
Temperature Up/Down. It also  
has a dual display which shows  
the SET Temp and the actual  
INSIDE Temp.

Automatic safety switch  
prevents overheating.

Fan provides circulation  
for even heat distribution.

Removable perforated  
bottom shelf allows for  
constant heat circulation  
as well as easy access to  
heating components.

This unit operates on  
110/120V AC and has the  
cord connected with a  
hospital grade cap.

FOR CART NUMBER	DESCRIPTION	SPACE BETWEEN SHELVES	OVERALL		
			L	W	H
110135	Stainless Steel Warming Cabinet	Stainless Steel Adjustable Wire Shelves	24.036"	30.125"	36.000"
108414-00	Stainless Steel Table with Crossbrace		36.000"	24.000"	34.000"



**Our Classic** single door warming cabinet is 36.00" tall by 30.125" wide and 24.036" deep. It is designed for the warming and storage of both post-surgery blankets and irrigation or intravenous solutions.

The compact size makes it perfect for areas with limited space and can be placed on an existing countertop, inserted into a recessed area or ordered with a mobile cart.

The easy to use, digital readout, allows for the user to set and maintain a specific temperature. An illuminated switch shows when

the unit is powered on, while additional lights and sounds alarm to prevent overheating. The cabinet is self contained with no heat vents. A fan circulates the warm air to maintain a constant temperature throughout the unit, making it more efficient and the controller more accurate. The temperature will be maintained with no heat leakage into the outside area, so the cabinet stays cool to the touch.

This unit has 6 positions available for wire shelves, allowing for custom separation of solutions and blankets.