

"LEADER IN HOSPITAL CARE PRODUCTS FOR OVER 60 YEARS"



**SSCI**

**SUBURBAN SURGICAL Co., Inc.**

275 Twelfth Street, Wheeling, IL 60090

All intellectual material, including the text, graphics and other components of this catalog are Copyrighted by SSCI.  
Any use or duplication of these intellectual materials, without written consent by SSCI, is strictly forbidden.

# RECEPTACLES

**PRECISION CRAFTED STAINLESS STEEL EQUIPMENT  
DESIGNED FOR LONG-LASTING PERFORMANCE**



## Materials:

### Chassis:

- Frame Members:** 16 gauge, 1" dia. tubular stainless steel, type 304, bright buffed finish
- Cross Members:** 16 gauge, 1" dia. tubular stainless steel, type 304, bright buffed finish
- Tray Supports:** .375" dia. rod stainless steel, type 304, bright buffed finish
- Casters:** 2" dia. non-marking, swivel ball bearing, zinc plated

## Catch Basin

Our stainless steel Catch Basin Receptacle is designed for convenience and ease of mobility. It includes a removable 15 quart Rectangular Tray.

P/N	DESCRIPTION
108130-00	SS Mobile Catch Basin, 15 Quart Rectangular Tray 26.810"L x 13.935"W x 20"H Overall
021950P	Additional Bin

## Construction:

### Chassis:

- Frame Members:** one piece "U"-shaped; bends are mandrill-bent for wrinkle-free bends
- Cross Members:** ends are coped for uniform fitting joints; fully heli-arc circumference-welded to Frame Members; welds are buffed to a smooth, bright polished finish
- Tray Supports:** fully heli-arc circumference-welded to Frame Members; welds are buffed to a smooth, bright polished finish
- Casters:** stem-mounted into caster socket with caster socket inserted into ends of Frame Member



**Mobile Sponge Receptacle with Cross Brace**



**Mobile Sponge Receptacle with Shelf**

## Materials:

### Chassis with Cross Brace:

- Frame Members:** 16 gauge, 1" dia. tubular stainless steel, type 304, bright buffed finish
- Pail Support:** 11 gauge, stainless steel, type 304, #4 finish
- Cross Braces:** .375" dia. rod stainless steel, type 304, bright buffed finish
- Connecting Rail:** .375" dia. rod stainless steel, type 304, bright buffed finish

### Chassis with Shelf:

- Frame Members:** 16 gauge, 1" dia. tubular stainless steel, type 304, bright buffed finish
- Pail Support:** 11 gauge, stainless steel, type 304, #4 finish
- Shelf:** 18 gauge, stainless steel, type 304, #4 finish
- Casters:** 3" dia. non-marking, swivel ball bearing, zinc plated

## Mobile Sponge Receptables

Our Mobile Sponge Receptacles are designed for convenience and ease of mobility and include a removable 13 quart Pail. They are available with either a sturdy Cross Brace or handy Shelf.

P/N	DESCRIPTION/DIMENSIONS
108152-00	SS Mobile Sponge Receptacle with Cross Brace 13 Quart Pail 14.130"L x 12.866"W x 34"H Overall
108153-00	SS Mobile Sponge Receptacle with Shelf 13 Quart Pail 14.130"L x 12.866"W x 34"H Overall
021876P	Additional 13 quart (12.30-liter) Bucket

## Construction:

### Chassis with Cross Brace or Shelf:

- Frame Members:** one piece "U"-shaped; bends are mandrill-bent for wrinkle-free bends
- Pail Support:** one piece formed ring; fully heli-arc welded to Frame Members; welds are buffed to a smooth, bright polished finish
- Cross Braces:** formed; ends fully heli-arc circumference-welded to Frame Members; welds are buffed to a smooth, bright polished finish
- Connecting Rail:** ends are fully heli-arc circumference welded to Cross Braces; welds are buffed to a smooth, bright polished finish
- Shelf:** one piece with sides hemmed; hemmed edges are flanged down creating a unitized shelf and double-thick angle frame; corners are notched to fit the Frame Member contour; unitized shelf and double-thick angle frame's vertical flanges are fully heli-arc welded to Frame Members (welds are concealed)
- Casters:** stem-mounted into caster socket with caster socket inserted into ends of Frame Member



## Mobile Sponge Receptacle

Stainless steel base features a durable heli-arc welded frame with molded rubber pads and four 2-in. (5.08cm) diameter non-marking, swivel casters. Securely holds basin to keep it close to your work area. Vinyl wrap-around bumper protects frame and facility walls from damage. Available in a 7-quart (6.62-liter) round basin or a 10-quart (9.46-liter) oval basin.

P/N	DESCRIPTION/DIMENSIONS
108150-00	Base with 7-quart (6.62-liter) Round Basin 17.12W X 6.98H X 17.12D in. (43.48W x 17.73H X 43.48D cm)
021758P	Additional 7-quart (6.62-liter) Round Basin
108145-00	Base with 10-quart (9.46-liter) Oval Basin 20.37W X 6.98H X 16.62D in. (51.74W x 17.73H X 42.21D cm)
021714V	Additional 10-quart (9.46-liter) Oval Basin



**Mobile Sponge Receptacle with Round Basin**



**Mobile Sponge Receptacle with Oval Basin**



## Mobile Kick Bucket

Durable heli-arc welded stainless steel frame with molded rubber pads holds a 13-quart (12.30-liter) stainless steel bucket securely. Non-marking 2-in. (5.08cm) diameter casters swivel for superior mobility. Vinyl wrap-around bumper protects frame and your facility's furniture and walls from damage. Measures 13.50-in. diameter x 14-in. high (34.29cm diameter x 35.56cm high) overall.

P/N	DESCRIPTION
108140-00	Includes 13-quart (12.30-liter) Bucket

## Kick Bucket

An indispensable component of any surgical suite. Heavy-duty stainless steel bucket may be used under collection holes on operating tables or as an extra or replacement bucket for Solution Stands and Kick Bucket. Swinging handle folds out of the way.

P/N	DESCRIPTION
021876P	13 quart (12.30-liter) Bucket



## Materials:

### Base:

- Frame Member:** 11 gauge, stainless steel, type 304, 2B finish
- Pail/Basin Supports:** 11 gauge, stainless steel, type 304, 2B finish
- Pail/Basin Pads:** molded rubber
- Cross Members:** 11 gauge, stainless steel, type 304, 2B finish
- Wraparound Bumper:** non-marring extruded vinyl
- Casters:** 2" dia. non-marking, swivel ball bearing, zinc plated

### Fasteners:

- Acorn Nuts:** stainless steel, type 18-8
- Cap Screws:** stainless steel, type 18-8

## Construction:

### Base:

- Frame Member:** one piece formed into a ring; ends are butted and heliarc welded (welds are concealed)
- Pail/Basin Supports:** four formed uprights punched and fitted with Pail/Basin Pads; punched to accommodate Casters; fully heli-arc welded to Frame Member; welds are polished to a smooth, bright polished finish
- Cross Members:** full-length; punched and fitted with Pail/Basin Pads; punched to accommodate Casters
- Wraparound Bumper:** continuous one piece; wrapped around Frame Member face and edges
- Casters:** threaded-stem mounted through holes in Cross Members and Pail/Basin Supports and secured with acorn nut