

“LEADER IN ANIMAL CARE PRODUCTS FOR OVER 60 YEARS”



SSCI

SUBURBAN SURGICAL Co., INC.

275 Twelfth Street, Wheeling, IL 60090

All intellectual material, including the text, graphics and other components of this catalog are Copyrighted by SSCI. Any use or duplication of these intellectual materials, without written consent by SSCI, is strictly forbidden.

Catalog No. 702488-0308

For Product Information Call TOLL-FREE 1-800-323-7366 • (847) 537-9320 • Fax (847) 537-9061
Hours: 8:30am - 5:00pm CST Monday - Friday
Visit our website at: www.suburbansurgical.com

ANCILLARY TABLES

**PRECISION CRAFTED STAINLESS STEEL EQUIPMENT
DESIGNED FOR LONG-LASTING PERFORMANCE**

Now you can outfit your entire facility with the best stainless steel Ancillary Tables from SSCI. The high-quality tables on the following pages are designed to be practical, functional and long lasting. Our goal is to offer you

the items you need to enhance the productivity of your facility.

Whether you're looking for Instrument Tables, Over Instrument Tables or Utility Tables, you'll find it at SSCI.



Instrument Table with Cross Brace



Over Instrument Table



Utility Tables



Instrument Table with Shelf

These Utility Tables feature built-in drawers to keep your instruments and treatment supplies clean, organized and easily accessible wherever you need them within your facility. Type 304 stainless steel construction offers long-lasting durability with easy care. Drawers are a full 4 in. (10.16cm) deep and feature a friction latch to stay closed while Table is being moved. Ball-bearing slides allow drawer to open fully, yet feature a positive stop to prevent drawers from being removed accidentally.

All welds are ground smooth and buffed to a smooth, bright polished finish.

Retaining rail around three sides of top surface help keep supplies in place during transit.

Tubular stainless steel legs are welded to drawer housing and shelves for outstanding strength.

Non-marking 3-in. (7.62cm) diameter ball-bearing swivel casters let you easily roll Table where it is needed.



1 Drawer Utility Table

P/N	OVERALL SIZE	
	Inches	Centimeters
108495-00	20.750W x 35.656H x 16.375D	52.705W x 90.57H x 41.59D

Instrument Tables

Choose our stainless steel Instrument Tables to keep your instruments organized and right at hand during procedures. Use it to lay out instruments in a central supply area, then easily roll it to where instruments will be used.

Constructed of Type 304 stainless steel for exceptionally easy care and long life. Choose from two styles; with a bottom cross brace or a shelf. Both styles are available in a choice of lengths and either with or without a drawer mounted beneath the top to offer greater storage capacity.

Tubular stainless steel legs are welded to top for outstanding strength.

Tubular stainless steel cross brace gives these tables outstanding rigidity while allowing for leg room underneath.

Non-marking 3-in. (7.62cm) diameter ball-bearing swivel casters let you easily roll Table where it is needed.

All welds are ground smooth and buffed to a smooth, bright polished finish.




Lower shelf lets you keep additional supplies right at hand.

Instrument Table

With Cross Brace	With Shelf	OVERALL SIZE	
		Inches	Centimeters
108410-00	108475-00	20W x 34H x 16D	50.80W x 86.36H x 40.64D
108414-00	108479-00	36W x 34H x 24D	91.44W x 86.36H x 60.96D
108415-00	108480-00	48W x 34H x 20D	121.92W x 86.36H x 50.80D
108416-00	108481-00	48W x 34H x 24D	121.92W x 86.36H x 60.96D

Over Instrument Table



Over Instrument Table

Use this Table to keep surgical instruments, supplies, even monitoring equipment conveniently close at hand during procedures. Manufactured entirely of durable Type 304 stainless steel, our Over Instrument Table features a mobile telescoping base for versatility. Plus, our conveniently-designed crank operates from either side of the Table to raise top height from 39-in. to 62-in. (99.06cm to 157.48cm). Table top surface measures a full 40-in. long x 23-in. wide (101.60cm long x 58.42cm wide) to hold an assortment of supplies. And the raised back rail keeps equipment from sliding off accidentally. Rolls smoothly on 3-in. (7.62cm) diameter non-marking casters to allow you to move Table into position easily. Two casters have brakes for safety.



Over Instrument Table can be raised and lowered to adjust to the height of your operating table.

P/N	DIMENSIONS	
	Inches	Centimeters
108431-23	40W x 39 to 62H x 23D	101.60W x 99.06 to 157.48H x 58.42D

For Product Information Call TOLL-FREE 1-800-323-7366 • (847) 537-9320 • Fax (847) 537-9061
Hours: 8:30am - 5:00pm CST Monday - Friday
 Visit our website at: www.suburbansurgical.com