

“LEADER IN ANIMAL CARE PRODUCTS FOR OVER 60 YEARS”



**SSCI**

**SUBURBAN SURGICAL Co., INC.**

275 Twelfth Street, Wheeling, IL 60090

All intellectual material, including the text, graphics and other components of this catalog are Copyrighted by SSCI. Any use or duplication of these intellectual materials, without written consent by SSCI, is strictly forbidden.

**Catalog No. 702488-0308**

**For Product Information Call TOLL-FREE 1-800-323-7366 • (847) 537-9320 • Fax (847) 537-9061**

**Hours: 8:30am - 5:00pm CST Monday - Friday**

Visit our website at: [www.suburbansurgical.com](http://www.suburbansurgical.com)

# ANCILLARY EQUIPMENT

**PRECISION CRAFTED STAINLESS STEEL EQUIPMENT  
DESIGNED FOR LONG-LASTING PERFORMANCE**

**Now** you can outfit your entire facility with the best stainless steel Ancillary Equipment from SSCI. The high-quality tables on the following pages are designed to be practical, functional and long lasting. Our goal is to offer you the items you need to enhance the productivity of your facility.

Whether you're looking for Mayo Stands, IV Stands, Solution Stands or Sponge Receptacles, you'll find it at SSCI.



**Mobile IV Stand**



**Mayo Stand**



**Solution Stand**



**Kick Bucket**



**Sponge Receptacle**



**Bucket**

**Sturdy,** all stainless steel Mayo Stands feature a heavy, weighted base to ensure unflinching balance and extraordinary stability. Choose our Regular, two-caster model with two extended support feet, or our Emergency model with three casters to allow you to move the Mayo Stand more quickly to where it is needed.

Both styles feature heli-arc welded polished joints for superior

strength and long lasting stability. The Type 304 stainless steel tray top offers plenty of space for supplies and removes easily for cleaning and sterilizing. Smooth-rolling casters make it easy to move our Mayo Stands throughout your facility.

Both models are offered in two tray top sizes in Foot-operated or Hand-operated styles to suit your needs.

**Emergency, Foot-Operated Model**



**Regular, Hand-Operated Model**



**Foot-Operated Mayo Stands**

Simply lift the tray up to the desired height from 38-in. to 61-in. (96.52cm to 154.94cm). It automatically locks in place. Then just depress the foot pedal to lower top when needed.

P/N FOOT OPERATED	TRAY SIZE (WxD)		OVERALL SIZE (WxD)	
	Inches	Centimeters	Inches	Centimeters
<b>Regular</b>				
108226-00	19.125 x 12.625	48.58 x 32.07	19.50 x 17.50	49.53 x 44.45
108228-00	21.250 x 16.250	53.98 x 41.28	21.625 x 17.50	54.93 x 44.45
<b>Emergency</b>				
108231-00	19.125 x 12.625	48.58 x 32.07	19.50 x 17.50	49.53 x 44.45
108233-00	21.250 x 16.250	53.98 x 41.28	21.625 x 17.50	54.93 x 44.45
<b>Additional Mayo Tray</b>				
021905P	19.125 x 12.625	48.58 x 32.07	19.50 x 17.50	49.53 x 44.45
021906P	21.250 x 16.250	53.98 x 41.28	21.625 x 17.50	54.93 x 44.45

**Hand-Operated Mayo Stands**

Knob-operated collet positively locks top at the desired height from 39-in. to 62-in. (99.06cm 157.48cm) high. Simply loosen knob to raise or lower top, then tighten to hold securely in place.

P/N HAND OPERATED	TRAY SIZE (WxD)		OVERALL SIZE (WxD)	
	Inches	Centimeters	Inches	Centimeters
<b>Regular</b>				
108246-00	19.125 x 12.625	48.58 x 32.07	20.875 x 17.50	53.02 x 44.45
108248-00	21.250 x 16.250	53.98 x 41.28	23.000 x 17.50	58.42 x 44.45
<b>Emergency</b>				
108251-00	19.125 x 12.625	48.58 x 32.07	20.875 x 17.50	53.02 x 44.45
108253-00	21.250 x 16.250	53.98 x 41.28	23.000 x 17.50	58.42 x 44.45
<b>Additional Mayo Tray</b>				
021905P	19.125 x 12.625	48.58 x 32.07	19.50 x 17.50	49.53 x 44.45
021906P	21.250 x 16.250	53.98 x 41.28	21.625 x 17.50	54.93 x 44.45

**For Product Information Call TOLL-FREE 1-800-323-7366 • (847) 537-9320 • Fax (847) 537-9061**  
**Hours: 8:30am - 5:00pm CST Monday - Friday**  
 Visit our website at: [www.suburbansurgical.com](http://www.suburbansurgical.com)

# Solution Stands and Sponge Receptacles



**Mobile Solution Stand with One Basin and Cross Brace**



**Mobile Solution Stand with Two Basins and Shelf**

## Mobile Solution Stands

Convenient Mobile Solution Stands are designed to keep 7-quart(s) (6.62-liter) round basin close at hand. Sturdy frame is made of 1-in. (2.54cm) diameter Type 304 stainless steel tubing, heli-arc welded at all joints for extra strength. Welds are buffed and polished to a bright shine. Comes complete with one or two 7-quart (6.62-liter) round basin(s) and four 3-in. (7.62cm) diameter non-marking, swivel casters.

P/N	DESCRIPTION/DIMENSIONS
108270-00	Stand with Cross Brace; one 7-qt. (6.62-liter) Round Basin 13.88W X 34.00H X 13.88D in. (35.26W X 86.36H X 35.26D cm)
108271-00	Stand with Shelf; one 7-qt. (6.62-liter) Round Basin 13.88W X 34.00H X 13.88D in. (35.26W X 86.36H X 35.26D cm)
108272-00	Stand with Cross Brace; two 7-qt. (6.62-liter) Round Basin 29.88W X 34.00H X 13.88D in. (75.90W X 86.36H X 35.26D cm)
108273-00	Stand with Shelf; two 7-qt. (6.62-liter) Round Basin 29.88W X 34.00H X 13.88D in. (75.90W X 86.36H X 35.26D cm)
021758P	Additional 7-quart (6.62-liter) Round Basin

## Mobile Sponge Receptacle

Stainless steel base features a durable heli-arc welded frame with molded rubber pads and four 2-in. (5.08cm) diameter non-marking, swivel casters. Securely holds basin to keep it close to your work area. Vinyl wrap-around bumper protects frame and facility walls from damage. Available in a 7-quart (6.62-liter) round basin or a 10-quart (9.46-liter) oval basin.

P/N	DESCRIPTION/DIMENSIONS
108150-00	Base with 7-quart (6.62-liter) Round Basin 17.12W X 6.98H X 17.12D in. (43.48W x 17.73H X 43.48D cm)
021758P	Additional 7-quart (6.62-liter) Round Basin
108145-00	Base with 10-quart (9.46-liter) Oval Basin 20.37W X 6.98H X 16.62D in. (51.74W x 17.73H X 42.21D cm)
021714V	Additional 10-quart (9.46-liter) Oval Basin



**Mobile Sponge Receptacle with Round Basin**



**Mobile Sponge Receptacle with Oval Basin**



## Mobile Kick Bucket

Durable heli-arc welded stainless steel frame with molded rubber pads holds a 13-quart (12.30-liter) stainless steel bucket securely. Non-marking 2-in. (5.08cm) diameter casters swivel for superior mobility. Vinyl wrap-around bumper protects frame and your facility's furniture and walls from damage. Measures 13.50-in. diameter x 14-in. high (34.29cm diameter x 35.56cm high) overall.

P/N	DESCRIPTION
108140-00	Includes 13-quart (12.30-liter) Bucket

## Bucket

An indispensable component of any surgical suite. Heavy-duty stainless steel bucket may be used under collection holes on operating tables or as an extra or replacement bucket for Solution Stands and Kick Bucket. Swinging handle folds out of the way.


P/N	DESCRIPTION
021876P	13 quart (12.30-liter) Bucket



**Our** heavy Type 304 stainless steel IV Stands feature a uniquely-designed heavy, weighted base to ensure unfailing balance and extraordinary stability. Stainless steel hanger loops are balanced to hold IV bags securely without tipping. Strong, durable heli-arc welded joints are buffed and polished for a finished look. And our smooth-rolling 2-in. (5.08cm) diameter non-marking cast-

ers make it easy to move our IV Stands throughout your facility. Choose from models with two or four IV hangers. All are available in Foot-operated or Hand-operated styles to suit your needs. Measure 21-in. long x 21-in. wide (53.34cm long x 53.34cm wide) overall.

## Foot-Operated Models




### Foot-Operated IV Stands

Lift the upper portion of this IV Stand to raise the hanger height from 61-in. to 101-in. (154.94cm to 256.54) and automatically lock in place. Foot pedal releases lock to lower hangers when needed.

FOOT OPERATED	DESCRIPTION
108180-00	2-Hook
108180-01	4-Hook

## Hand-Operated Models



### Hand-Operated IV Stands

Knob-operated collet positively locks hangers at any desired height from 61-in. to 101-in. (154.94cm to 256.54) high. Simply loosen knob to change top height, then tighten to hold securely in place.

HAND OPERATED	DESCRIPTION
108200-00	2-Hook
108200-01	4-Hook